
Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

GK2

Uniprot

Q14410

Synonyms:

GK2; GKP2; GKTA; Glycerol kinase 2; GK 2; Glycerokinase 2; ATP: glycerol 3-phosphotransferase 2; Glycerol kinase; testis specific 2


Immunogen:

Synthesized peptide derived from the Internal region of human GK2.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 A549 cells using GK2 Polyclonal Antibody.



Western Blot analysis of Jurkat cells using GK2 Polyclonal Antibody.