

GPR152 Antibody



PACO06525

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

GPR152

Uniprot

Q8TDT2

Synonyms:

GPR152; PGR5; Probable G-protein coupled receptor 152; G-protein coupled receptor PGR5


Immunogen:

Synthesized peptide derived from the N-terminal region of human GPR152.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of various cells using GPR152 Polyclonal Antibody.



Western Blot analysis of k562 cells using GPR152 Polyclonal Antibody.