

HDAC3 Antibody



PACO06538

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

HDAC3

Uniprot

O15379

Synonyms:

HDAC3; Histone deacetylase 3; HD3; RPD3-2; SMAP45


Immunogen:

Synthesized peptide derived from human HDAC3 around the non-phosphorylation site of S424.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HELA cells using HDAC3 Polyclonal Antibody.



Western Blot analysis of hela cells using HDAC3 Polyclonal Antibody.