

# HDAC4 Antibody



PACO06539

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:20000, WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

HDAC4

**Uniprot**

P56524

**Synonyms:**

HDAC4; KIAA0288; Histone deacetylase 4; HD4

**Immunogen:**


Synthesized peptide derived from human HDAC4 around the non-phosphorylation site of S632.

**Storage:**


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of HELA cells using HDAC4 Polyclonal Antibody.



Western Blot analysis of HELA cells using HDAC4 Polyclonal Antibody.