

HDAC5 Antibody



PACO06540

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

HDAC5

Uniprot

Q9UQL6

Synonyms:

HDAC5; KIAA0600; Histone deacetylase 5; HD5; Antigen NY-CO-9


Immunogen:

Synthesized peptide derived from human HDAC5 around the non-phosphorylation site of S498.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 NIH-3T3 cells using HDAC5 Polyclonal Antibody.



Western Blot analysis of K562 cells using HDAC5 Polyclonal Antibody.