

HDAC5/HDAC9 Antibody



PACO06541

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

HDAC5/HDAC9

Uniprot

Q9UQL6/Q9UKV0

Synonyms:

HDAC5; KIAA0600; Histone deacetylase 5; HD5; Antigen NY-CO-9; HDAC9; HDAC7; HDAC7B; HDRP; KIAA0744; MITR; Histone deacetylase 9; HD9; Histone deacetylase 7B; HD7; HD7b; Histone deacetylase-related protein; MEF2-interacting transcription rep


Immunogen:

Synthesized peptide derived from human HDAC5/9 around the non-phosphorylation site of S259/220.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 RAT-intestines cells using HDAC5/9 Polyclonal Antibody.



Western Blot analysis of HepG2 cells using HDAC5/9 Polyclonal Antibody.