

# HDAC8 Antibody



PACO06543

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat, Monkey

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:20000, WB:1:500-1:2000,  
IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

HDAC8

**Uniprot**

Q9BY41

**Synonyms:**

HDAC8; HDACL1; CDA07; Histone deacetylase 8; HD8

**Immunogen:**


Synthesized peptide derived from human HDAC8 around the non-phosphorylation site of S39.

**Storage:**


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of NIH-3T3 HeLa Jurkat cells using HDAC8 Polyclonal Antibody.



Western Blot analysis of 293-UV cells using HDAC8 Polyclonal Antibody.