

# HDAC9 Antibody



PACO06548

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000,  
IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:****Gene ID:**

HDAC9

**Uniprot**

Q9UKV0-5,6,7

**Synonyms:**

Histone deacetylase 9; HD9; Histone deacetylase 7B; HD7; HD7b; Histone deacetylase-related protein; MEF2-interacting transcription repressor MITR; HDAC9; HDAC7; HDAC7B; HDRP; KIAA0744; MITR

**Immunogen:**


Synthesized peptide derived from the C-terminal region of human Histone deacetylase 9.

**Storage:**


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of K562 cells using Histone deacetylase 9 Polyclonal Antibody.



Western Blot analysis of HepG2 cells using Histone deacetylase 9 Polyclonal Antibody.