

# HNRNPA1 Antibody



PACO06556

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

ELISA:1:20000, WB:1:500-1:2000,  
IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:****Gene ID:**

HNRNPA1

**Uniprot**

P09651

**Synonyms:**

HNRNPA1; HNRPA1; Heterogeneous nuclear ribonucleoprotein A1; hnRNP A1; Helix-destabilizing protein; Single-strand RNA-binding protein; hnRNP core protein A1

**Immunogen:**


Synthesized peptide derived from the Internal region of human hnRNP A1.

**Storage:**


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of Jurkat NIH-3T3 cells using hnRNP A1 Polyclonal Antibody.



Western Blot analysis of COLO205 cells using hnRNP A1 Polyclonal Antibody.