

HNRNPF Antibody



PACO06557

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

HNRNPF

Uniprot

P52597

Synonyms:

HNRNPF; HNRPF; Heterogeneous nuclear ribonucleoprotein F; hnRNP F; Nucleolin-like protein mcs94-1


Immunogen:

Synthesized peptide derived from the N-terminal region of human hnRNP F.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of various cells using hnRNP F Polyclonal Antibody.



Western Blot analysis of HepG2 cells using hnRNP F Polyclonal Antibody.