

HOXA11/HOXD11 Antibody



PACO06560

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

HOXA11/HOXD11

Uniprot

P31270/P31277

Synonyms:

HOXA11; HOX11; Homeobox protein Hox-A11; Homeobox protein Hox-11; HOXD11; HOX4F; Homeobox protein Hox-D11; Homeobox protein Hox-4F


Immunogen:

Synthesized peptide derived from the C-terminal region of human HoxA11/D11.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of LOVO cells using HoxA11/D11 Polyclonal Antibody.



Western Blot analysis of 293T-UV cells using HoxA11/D11 Polyclonal Antibody.