

HOXB9 Antibody



PACO06561

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

HOXB9

Uniprot

P17482

Synonyms:

HOXB9; HOX2E; Homeobox protein Hox-B9; Homeobox protein Hox-2.5; Homeobox protein Hox-2E


Immunogen:

Synthesized peptide derived from the C-terminal region of human Hox-B9.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HELA cells using Hox-B9 Polyclonal Antibody.



Western Blot analysis of HeLa cells using Hox-B9 Polyclonal Antibody.