

HSPA9 Antibody



PACO06566

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

HSPA9

Uniprot

P38646/Q8N1C8

Synonyms:

HSPA9; GRP75; HSPA9B; Stress-70 protein; mitochondrial; 75 kDa glucose-regulated protein; GRP-75; Heat shock 70 kDa protein 9; Mortalin; MOT; Peptide-binding protein 74; PBP74


Immunogen:

Synthesized peptide derived from the C-terminal region of human HSP A9.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 RAT-MUSCLE cells using HSP A9 Polyclonal Antibody.



Western Blot analysis of COS7 cells using HSP A9 Polyclonal Antibody.