

DLGAP5 Antibody



PACO06570

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:**Gene ID:**

DLGAP5

Uniprot

Q15398

Synonyms:

DLGAP5; DLG7; KIAA0008; Disks large-associated protein 5; DAP-5; Discs large homolog 7; Disks large-associated protein DLG7; Hepatoma up-regulated protein; HURP


Immunogen:

Synthesized peptide derived from the C-terminal region of human HURP.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 NIH-3T3 cells using HURP Polyclonal Antibody.



Western Blot analysis of A549 cells using HURP Polyclonal Antibody.