

IL17RC Antibody



PACO06582

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:

Gene ID:

IL17RC

Uniprot

Q8NAC3

Synonyms:

IL17RC; Interleukin-17 receptor C; IL-17 receptor C; IL-17RC; Interleukin-17 receptor homolog; IL17Rhom; Interleukin-17 receptor-like protein; IL-17RL; ZcytoR14


Immunogen:

Synthesized peptide derived from the C-terminal region of human IL-17RC.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of LOVO cells using IL-17RC Polyclonal Antibody.



Western Blot analysis of HuvEc cells using IL-17RC Polyclonal Antibody.