

# JAK2 Antibody



PACO06595

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000,  
IHC:1:100-1:300

**Protein Background:****Gene ID:**

JAK2

**Uniprot**

O60674

**Synonyms:**

JAK2; Tyrosine-protein kinase JAK2; Janus kinase 2; JAK-2

**Immunogen:**


Synthesized peptide derived from human JAK2 around the non-phosphorylation site of Y1007.

**Storage:**


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of PC3 KB NIH-3T3 cells using JAK2 Polyclonal Antibody.



Western Blot analysis of HT29 cells using JAK2 Polyclonal Antibody.