

MCL1 Antibody



PACO06628

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

MCL1

Uniprot

Q07820

Synonyms:

MCL1; BCL2L3; Induced myeloid leukemia cell differentiation protein Mcl-1; Bcl-2-like protein 3; Bcl2-L-3; Bcl-2-related protein EAT/mcl1; mcl1/EAT


Immunogen:

Synthesized peptide derived from the Internal region of human Mcl-1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 Jurkat cells using Mcl-1 Polyclonal Antibody.



Western Blot analysis of Jurkat cells using Mcl-1 Polyclonal Antibody.