

MTNR1B Antibody



PACO06633

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000

Protein Background:

Gene ID:

MTNR1B

Uniprot

P49286

Synonyms:

MTNR1B; Melatonin receptor type 1B; Mel-1B-R; Mel1b receptor


Immunogen:

Synthesized peptide derived from the C-terminal region of human MEL-1B-R.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HELA A549 cells using MEL-1B-R Polyclonal Antibody.



Western Blot analysis of HepG2 cells using MEL-1B-R Polyclonal Antibody.