

# MERTK/TYRO3 Antibody



PACO06637

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000

**Protein Background:****Gene ID:**

MERTK/TYRO3

**Uniprot**

Q12866/Q06418

**Synonyms:**

MERTK; MER; Tyrosine-protein kinase Mer; Proto-oncogene c-Mer; Receptor tyrosine kinase MerTK; TYRO3; BYK; DTK; RSE; SKY; Tyrosine-protein kinase receptor TYRO3; Tyrosine-protein kinase DTK; Tyrosine-protein kinase RSE; Tyrosine-protein kin

**Immunogen:**


Synthesized peptide derived from human MerTK/Tyro3 around the non-phosphorylation site of Y753/Y685.

**Storage:**


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of HuvEc cells using MerTK/Tyro3 Polyclonal Antibody.



Western Blot analysis of HuvEc cells using MerTK/Tyro3 Polyclonal Antibody.