

GRM4 Antibody



PACO06638

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

GRM4

Uniprot

Q14833

Synonyms:

GRM4; GPRC1D; MGLUR4; Metabotropic glutamate receptor 4; mGluR4


Immunogen:

Synthesized peptide derived from the C-terminal region of human mGluR-4.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAT-KIDNEY cells using mGluR-4 Polyclonal Antibody.



Western Blot analysis of 293 cells using mGluR-4 Polyclonal Antibody.