

MGMT Antibody



PACO06641

Product Information

Size:

50ug

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

MGMT

Uniprot

P16455

Synonyms:

MGMT; Methylated-DNA--protein-cysteine methyltransferase; 6-O-methylguanine-DNA methyltransferase; MGMT; O-6-methylguanine-DNA-alkyltransferase


Immunogen:

Synthesized peptide derived from the N-terminal region of human MGMT.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of MCF7 PC3 cells using MGMT Polyclonal Antibody.



Western Blot analysis of Jurkat cells using MGMT Polyclonal Antibody.