

MKNK1 Antibody



PACO06647

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

MKNK1

Uniprot

Q9BUB5

Synonyms:

MKNK1; MNK1; MAP kinase-interacting serine/threonine-protein kinase 1; MAP kinase signal-integrating kinase 1; MAPK signal-integrating kinase 1; Mnk1


Immunogen:

Synthesized peptide derived from the Internal region of human Mnk1.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of MOUSE-BRAIN cells using Mnk1 Polyclonal Antibody.



Western Blot analysis of HuvEc cells using Mnk1 Polyclonal Antibody.