

MEOX2 Antibody



PACO06649

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

MEOX2

Uniprot

P50222

Synonyms:

MEOX2; GAX; MOX2; Homeobox protein MOX-2; Growth arrest-specific homeobox; Mesenchyme homeobox 2


Immunogen:

Synthesized peptide derived from the Internal region of human MOX-2.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAT-MUSCLE cells using MOX-2 Polyclonal Antibody.



Western Blot analysis of COLO205 cells using MOX-2 Polyclonal Antibody.