

MRPS34 Antibody



PACO06657

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

MRPS34

Uniprot

P82930

Synonyms:

MRPS34; 28S ribosomal protein S34; mitochondrial; MRP-S34; S34mt


Immunogen:

Synthesized peptide derived from the C-terminal region of human MRP-S34.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat HuvEc cells using MRP-S34 Polyclonal Antibody.



Western Blot analysis of Jurkat cells using MRP-S34 Polyclonal Antibody.