

MYD88 Antibody



PACO06662

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300, IHC:1:200-1:1000

Protein Background:**Gene ID:**

MYD88

Uniprot

Q99836

Synonyms:

MYD88; Myeloid differentiation primary response protein MyD88


Immunogen:

Synthesized peptide derived from the Internal region of human MyD88.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAT-MUSCLE cells using MyD88 Polyclonal Antibody.



Western Blot analysis of COLO205 cells using MyD88 Polyclonal Antibody.