

MT-ND1 Antibody



PACO06673

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

MT-ND1

Uniprot

P03886

Synonyms:

MT-ND1; MTND1; NADH1; ND1; NADH-ubiquinone oxidoreductase chain 1; NADH dehydrogenase subunit 1


Immunogen:

Synthesized peptide derived from the Internal region of human ND1.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of L929 cells using ND1 Polyclonal Antibody.



Western Blot analysis of COLO205 cells using ND1 Polyclonal Antibody.