

NDUFB9 Antibody



PACO06674

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

NDUFB9

Uniprot

Q9Y6M9

Synonyms:

NDUFB9; LYRM3; UQOR22; NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 9; Complex I-B22; CI-B22; LYR motif-containing protein 3; NADH-ubiquinone oxidoreductase B22 subunit


Immunogen:

Synthesized peptide derived from the Internal region of human NDUFB9.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HELA cells using NDUFB9 Polyclonal Antibody.



Western Blot analysis of 293 cells using NDUFB9 Polyclonal Antibody.