

PACO06689

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

NGF

Uniprot

P01138

Synonyms:

NGF; NGFB; beta-nerve growth factor; beta-NGF


Immunogen:

Synthesized peptide derived from the Internal region of human NGF.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of various cells using NGF Polyclonal Antibody.



Western Blot analysis of Jurkat cells using NGF Polyclonal Antibody.