

NEMF Antibody



PACO06703

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

NEMF

Uniprot

O60524

Synonyms:

NEMF; SDCCAG1; Nuclear export mediator factor NEMF; Antigen NY-CO-1;
Serologically defined colon cancer antigen 1


Immunogen:

Synthesized peptide derived from the C-terminal region of human NY-CO-1.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 cells using NY-CO-1 Polyclonal Antibody.



Western Blot analysis of 3T3 cells using NY-CO-1 Polyclonal Antibody.