

OR10AG1 Antibody



PACO06708

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR10AG1

Uniprot

Q8NH19

Synonyms:

OR10AG1; Olfactory receptor 10AG1; Olfactory receptor OR11-160


Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 10AG1.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat K562 cells using Olfactory receptor 10AG1 Polyclonal Antibody.



Western Blot analysis of Jurkat cells using Olfactory receptor 10AG1 Polyclonal Antibody.