

CDKN1C Antibody



PACO06730

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

CDKN1C

Uniprot

P49918

Synonyms:

CDKN1C; KIP2; Cyclin-dependent kinase inhibitor 1C; Cyclin-dependent kinase inhibitor p57; p57Kip2


Immunogen:

Synthesized peptide derived from the N-terminal region of human p57.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAT-MUSCLE cells using p57 Polyclonal Antibody.



Western Blot analysis of rat-muscle cells using p57 Polyclonal Antibody.