

PDK1 Antibody



PACO06740

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000

Protein Background:**Gene ID:**

PDK1

Uniprot

Q15118

Synonyms:

PDK1; PDHK1; [Pyruvate dehydrogenase [lipoamide]] kinase isozyme 1; mitochondrial; Pyruvate dehydrogenase kinase isoform 1


Immunogen:

Synthesized peptide derived from the Internal region of human PDK1.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc cells using PDK1 Polyclonal Antibody.



Western Blot analysis of HeLa cells using PDK1 Polyclonal Antibody.