

PDPK1 Antibody



PACO06741

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

PDPK1

Uniprot

O15530

Synonyms:

PDPK1; PDK1; 3-phosphoinositide-dependent protein kinase 1; hPDK1


Immunogen:

Synthesized peptide derived from human PDK1 around the non-phosphorylation site of S241.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HELA COLO 293T cells using PDK1 Polyclonal Antibody.



Western Blot analysis of COLO205 cells using PDK1 Polyclonal Antibody.