

PIAS2 Antibody



PACO06751

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA: 1:20000

WB: 1:500-1:2000

Form:

Liquid

Product Type:

Polyclonal Antibody

Conjugate:

Non-conjugated

Gene ID:

PIAS2

Uniprot

O75928

Synonyms:

PIAS2; PIASX; E3 SUMO-protein ligase PIAS2; Androgen receptor-interacting protein 3; ARIP3; DAB2-interacting protein; DIP; Msx-interacting zinc finger protein; Miz1; PIAS-NY protein; Protein inhibitor of activated STAT x; Protein inhibitor

Immunogen:

Synthesized peptide derived from the N-terminal region of human PIASx.

Immunogen Species:

Human

Storage Buffer:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.


Storage:

Upon receipt, store at -20 °C or -80 °C. Avoid repeated freeze


Purification Method:

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Product Images



Western Blot analysis of NIH-3T3 Jurkat K562 cells using PIASx Polyclonal Antibody.



Western Blot analysis of COS7 cells using PIASx Polyclonal Antibody.