

PRKACA/PRKACB/PRKACG Antibody



PACO06760

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

PRKACA/PRKACB/PRKACG

Uniprot

P17612/P22694/P22612

Synonyms:

PRKACA; PKACA; cAMP-dependent protein kinase catalytic subunit alpha; PKA C-alpha; PRKACB; cAMP-dependent protein kinase catalytic subunit beta; PKA C-beta; PRKACG; cAMP-dependent protein kinase catalytic subunit gamma; PKA C-gamma


Immunogen:

Synthesized peptide derived from human PKA α ; / β ; / γ ; cat around the non-phosphorylation site of T198.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 22RV1 HELA cells using PKA alpha / beta / gamma cat Polyclonal Antibody.



Western Blot analysis of COLO205 cells using PKAalpha / beta / gamma cat Polyclonal Antibody.