

PROKR1 Antibody



PACO06765

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

PROKR1

Uniprot

Q8TCW9

Synonyms:

PROKR1; GPR73; PKR1; Prokineticin receptor 1; PK-R1; G-protein coupled receptor 73; G-protein coupled receptor ZAQ; GPR73a


Immunogen:

Synthesized peptide derived from the N-terminal region of human PKR1.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using PKR1 Polyclonal Antibody.



Western Blot analysis of COLO205 cells using PKR1 Polyclonal Antibody.