

PPARG Antibody



PACO06774

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

PPARG

Uniprot

P37231

Synonyms:

PPARG; NR1C3; Peroxisome proliferator-activated receptor gamma; PPAR-gamma;
Nuclear receptor subfamily 1 group C member 3


Immunogen:

Synthesized peptide derived from human PPAR- γ around the non-phosphorylation site of S112.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of various cells using PPAR- γ Polyclonal Antibody.



Western Blot analysis of HuvEc cells using PPAR- γ Polyclonal Antibody.