

RAB20 Antibody



PACO06784

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

RAB20

Uniprot

Q9NX57

Synonyms:

RAB20; Ras-related protein Rab-20


Immunogen:

Synthesized peptide derived from the Internal region of human Rab 20.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat K562 NIH-3T3 cells using Rab 20 Polyclonal Antibody.



Western Blot analysis of A549 cells using Rab 20 Polyclonal Antibody.