

RPL17 Antibody



PACO06797

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

RPL17

Uniprot

P18621

Synonyms:

RPL17; 60S ribosomal protein L17; 60S ribosomal protein L23; PD-1


Immunogen:

Synthesized peptide derived from the Internal region of human Ribosomal Protein L17.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of customer's (cat sample) using Ribosomal Protein L17 Polyclonal Antibody.



Western Blot analysis of COS7 cells using Ribosomal Protein L17 Polyclonal Antibody.