

RNF113B Antibody



PACO06803

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

RNF113B

Uniprot

Q8IZP6

Synonyms:

RNF113B; RNF161; ZNF183L1; RING finger protein 113B; Zinc finger protein 183-like 1


Immunogen:

Synthesized peptide derived from the Internal region of human RNF113B.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using RNF113B Polyclonal Antibody.



Western Blot analysis of 293 cells using RNF113B Polyclonal Antibody.