

TOP2A Antibody



PACO06837

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

TOP2A

Uniprot

P11388

Synonyms:

TOP2A; TOP2; DNA topoisomerase 2-alpha; DNA topoisomerase II; alpha isozyme


Immunogen:

Synthesized peptide derived from the N-terminal region of human Topo II α .


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using Topo II alpha Polyclonal Antibody.



Western Blot analysis of COS7 cells using Topo IIalpha Polyclonal Antibody.