

TXNRD2 Antibody



PACO06843

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

TXNRD2

Uniprot

Q9NNW7

Synonyms:

TXNRD2; KIAA1652; TRXR2; Thioredoxin reductase 2; mitochondrial; Selenoprotein Z; SelZ; TR-beta; Thioredoxin reductase TR3


Immunogen:

Synthesized peptide derived from the C-terminal region of human TrxR2.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAT-MUSCLE cells using TrxR2 Polyclonal Antibody.



Western Blot analysis of HT29 cells using TrxR2 Polyclonal Antibody.