

VEGFB Antibody



PACO06852

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

VEGFB

Uniprot

P49765

Synonyms:

VEGFB; VRF; Vascular endothelial growth factor B; VEGF-B; VEGF-related factor; VRF


Immunogen:

Synthesized peptide derived from the Internal region of human VEGF-B.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO 293T HELA cells using VEGF-B Polyclonal Antibody.



Western Blot analysis of Jurkat cells using VEGF-B Polyclonal Antibody.