

CD55 Antibody



PACO06870

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:

Gene ID:

CD55

Uniprot

P08174

Synonyms:

CD55; CR; DAF; Complement decay-accelerating factor; CD55


Immunogen:

Synthesized peptide derived from the Internal region of human CD55.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using CD55 Polyclonal Antibody.