

ALPP/ALPPL2 Antibody



PACO06876

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

ALPP/ALPPL2

Uniprot

P05187/P10696

Synonyms:

ALPP; PLAP; Alkaline phosphatase, placental type; Alkaline phosphatase Regan isozyme; Placental alkaline phosphatase 1; PLAP-1; ALPPL2; ALPPL; Alkaline phosphatase, placental-like


Immunogen:

Synthesized peptide derived from the N-terminal region of human ALPP/ALPPL2.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of customer's sample using ALPP/ALPPL2 Polyclonal Antibody.



Western Blot analysis of 293 cells using ALPP/ALPPL2 Polyclonal Antibody.