

**TUBA1A/TUBA1B/TUBA1C/TUBA3C/TUBA4A
Antibody
PACO06879**



Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:

Gene ID:

TUBA1A/TUBA1B/TUBA1C/TUBA3C/TUBA4A

Uniprot

Q71U36/P68363/Q9BQE3/Q13748/P68366

Synonyms:

TUBA1A; TUBA3; Tubulin alpha-1A chain; alpha-tubulin 3; Tubulin B-alpha-1; Tubulin alpha-3 chain; TUBA1B; Tubulin alpha-1B chain; alpha-tubulin ubiquitous


Immunogen:

Synthesized peptide derived from human Tubulin α around the non-acetylation site of K163.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562, PC12 cells using Tubulin alpha Polyclonal Antibody.