

# ITGA6 Antibody



PACO06888

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000

**Protein Background:****Gene ID:**

ITGA6

**Uniprot**

P23229

**Synonyms:**

ITGA6; Integrin alpha-6; CD49 antigen-like family member F; VLA-6; CD49f

**Immunogen:**


Synthesized peptide derived from the Internal region of human Integrin &alpha; 6.

**Storage:**


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of NIH-3T3 cells using Integrin alpha 6 Polyclonal Antibody.



Western Blot analysis of 3T3 cells using Integrin alpha 6 Polyclonal Antibody.