

PVRL1 Antibody



PACO06896

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

NECTIN1

Uniprot

Q15223

Synonyms:

PVRL1; HVEC; PRR1; Poliovirus receptor-related protein 1; Herpes virus entry mediator C; Herpesvirus entry mediator C; HveC; Herpesvirus Ig-like receptor; HlgR; Nectin-1; CD111


Immunogen:

Synthesized peptide derived from the Internal region of human Nectin 1.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of KB cells using Nectin 1 Polyclonal Antibody.



Western Blot analysis of KB cells using Nectin 1 Polyclonal Antibody.