

IL13RA2 Antibody



PACO06902

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:

Gene ID:

IL13RA2

Uniprot

Q14627

Synonyms:

IL13RA2; IL13R; Interleukin-13 receptor subunit alpha-2; IL-13 receptor subunit alpha-2; IL-13R subunit alpha-2; IL-13R-alpha-2; IL-13RA2; Interleukin-13-binding protein; CD213a2


Immunogen:

Synthesized peptide derived from the Internal region of human IL-13R α 2.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2(UV) L929 cells using IL-13R alpha 2 Polyclonal Antibody.



Western Blot analysis of HEPG2-UV cells using IL-13Ralpha 2 Polyclonal Antibody.