

LEFTY1/LEFTY2 Antibody



PACO06923

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

LEFTY1/LEFTY2

Uniprot

O75610/O00292

Synonyms:

LEFTY1; LEFTB; LEFTYB; Left-right determination factor 1; Left-right determination factor B; Protein lefty-1; Protein lefty-B; LEFTY2; EBAF; LEFTA; LEFTYA; TGFB4; PSEC0024; Left-right determination factor 2; Endometrial bleeding-associated factor; Left-ri


Immunogen:

Synthesized peptide derived from the C-terminal region of human Lefty.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2-UV cells using Lefty Polyclonal Antibody.



Western Blot analysis of HEPG2-UV cells using Lefty Polyclonal Antibody.